Journée Alpview 14 novembre 2024

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
API de Xact

- Each object can be controled through attributes, commands and requests:
 - The **attributes** correspond to the **parameters** of the machine
 - Example : `http://localhost/api/Imager/attributes/avrFrames`
 - The **commands** correspond to **actions** on the machine
 - Example : `http://localhost/api/Imager/command/start`
 - The **requests** permit:
 - to retrieve the status and errors of each object
 - some special operations to start a live imager broadcasting or to get reconstructed data

API operations	HTTP		Websocket
	GET	POST	
Commands		Execute an action on the machine	
Attributes	Return the value of the parameter	Change the value of the parameter	
Requests	Return the status and the errors of the machine or object		Imager streaming / UniCT data retriever

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Permission Token

2.4. Permission token to modify the resources

Since concurrent access to modify the machine may happen and lead to hazardous behaviour, two special commands have to be used to control the machine, **hold** and **release**. The first command requests a permission token that allows the client to change the state of the machine. The second command releases the token once the client has finished to modify the machine state.

Note: only one token can be provided at the same time. The python script below shows an example of how to get and release such token.

Hold a permission token

Python script

```
import requests
import json

response = requests.get("https://XACT_SERVER_ADDRESS/API/hold")
json_response = json.loads(response.text)
```

Content of the json_response variable

```
{
  "object": "",
  "operation": "hold",
  "response": {
    "type": "message",
    "message": "07d438945d16459daa7b22be7dd8e23b"
  }
}
```

Release a permission token

Python script

```
import requests
import json

post_body = {"permission": "07d438945d16459daa7b22be7dd8e23b", "parameters": []}
response = requests.post("https://XACT_SERVER_ADDRESS/API/release", json = post_body)
json_response = json.loads(response.text)
```

Content of the json_response variable

```
{
  "object": "",
  "operation": "release",
  "response": {
    "type": "message",
    "message": "The machine has been successfully released"
  }
}
```

[Home](#)[API errors](#)[API](#)[Auto](#)[Automation](#)[AxisControl](#)[Console](#)[DisplayZoom](#)[Equalizer](#)[GeometryExporter](#)[GeometrySolver](#)[Guided](#)[Horizontal](#)[ImageCalculator](#)[Imager](#)[MacroEditor](#)[MainWindowInterface](#)[Remote](#)[Rotation](#)[RunScript](#)

1. Introduction

This API is a web service that permits to access and program RX Solutions X-Act software to perform a well-defined sequence of actions. This web service proposes several ways to communicate with X-Act:

- GET method to get the current state of the machine (cf. sections 2 and 3),
- POST method to change the state of the machine (cf. sections 2 and 3),
- WebSocket protocol to stream imager live or to receive reconstructed data (cf. section 4).

There are several resources available on the machine (xRay sources, imagers, axes, ...) listed on the left menu. For each resource a detailed list of the available attributes, commands and requests is provided.

2. Special commands

2.1. List all available resources and all available attributes, commands and requests of a resource

It is possible to have a list of all available resources that can be controlled through the API, and, for each resource, it is also possible to obtain a list of all attributes, commands and requests. The python script below shows an example of how to get these lists.

Get all available resources

Python script

```
import requests
import json

response = requests.get("https://XACT_SERVER_ADDRESS/API")
json_response = json.loads(response.text)
```

Content of the json_response variable

```
{
  "object": "",
  "operation": "requestAvailableObjects",
  "response": {
    "message": {
      "Axis": ["Horizontal", "ImagerH", "ImagerV", "ImagerZoom", "Rotation", "Vertical", "Zoom"],
      "Imager": ["Imager", "Lucas"],
      "Io": [{"Io": "Io"}],
      "Misc": ["Tomography", "INach3", "KeepOut", "Laminography", "Automation", "BeckhoffTemperatureSensor", "DisplayZoom", "Ge"],
      "Source": [{"Xray150"}],
      "TouchLcd": [{"ILcd"}]
    },
    "type": "message"
  }
}
```

Get available attributes, commands and requests for Xray150 resource

Python script

```
import requests
import json

response = requests.get("https://XACT_SERVER_ADDRESS/API/Xray150")
json_response = json.loads(response.text)
```

Content of the json_response variable

```
{
  "object": "Xray150",
  "operation": "requestObjectInfos",
  "response": {
    "message": {
      "attributes": [{"name": "connected", "type": "bool"}, {"name": "current", "type": "int"}, {"name": "filter", "type": "QGe"},
      "commands": [{"name": "requestSpotsize", "parametersType": "double", "returnType": "bool"}, {"name": "setCurrent", "param"},
      "states": [{"state": "HM_NOT_READY", "value": 1}, {"state": "VOLTAGE", "value": 2}, {"state": "CURRENT", "value": 4}, {"s"},
      "type": "message"
    }
  }
}
```

2.2. Error log

A global error log is available and permits to trace the event errors on the machine. This log can also be called for each resource as a request (cf. *Requests* section of each resource).

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

Page [Discussion](#)

Design pattern

A **design pattern**, also know as a software design pattern, is a reusable solution to a software eng templates that you and others can reuse for future projects.

Described below some useful design patterns that a developer can use in their application architect

Basic Design Patterns

State Machine

See more on the [State Machine Design Pattern](#).

Event Handler

See more on the [Event Handler Design Pattern](#).

Master/Slave

See more on the [Master/Slave Design Pattern](#).

Producer/Consumer

See more on the [Producer/Consumer Design Pattern](#).

Intermediate Design Patterns

Queued Message Handler (QMH)

See more on the [Queued Message Handler \(QMH\) Design Pattern](#).

Queued State Machine (QSM)

See more on the [Queued State Machine \(QSH\) Design Pattern](#).

Action Engine (AE)

See more on the [Action Engine \(AE\) Design Pattern](#) a.k.a. [Functional global variable](#).

Advanced Design Patterns

Object-Oriented Design Patterns

See more on [Object-Oriented Design Patterns](#).

Actor Oriented Design

See more on the [Actor Oriented Design Patterns](#).