

Bremsstrahlung and Ion Beam Current Measurements with SuSI ECR Ion Source

Tommi Ropponen
NSCL/MSU
East Lansing, MI, USA
ropponen@nscl.msu.edu

August 25th 2010, ECRIS'10, Grenoble, France




OVERVIEW

1. Measurement setup


- SuSI
- Measurement geometry
- DAQ

2. Different collimation schemes




OVERVIEW

1. Measurement setup
 - SuSI
 - Measurement geometry
 - DAQ
2. Different collimation schemes
3. FlatB vs. Bmin operation




OVERVIEW

1. Measurement setup
 - SuSI
 - Measurement geometry
 - DAQ
2. Different collimation schemes
3. FlatB vs. Bmin operation
4. The effect of ∇B
 - With fixed plasma length and B_{min}




OVERVIEW

1. Measurement setup
 - SuSI
 - Measurement geometry
 - DAQ
2. Different collimation schemes
3. FlatB vs. Bmin operation
4. The effect of ∇B
 - With fixed plasma length and $B_{\min}$
5. Summary




OVERVIEW

1. Measurement setup
 - SuSI
 - Measurement geometry
 - DAQ
2. Different collimation schemes
3. FlatB vs. Bmin operation
4. The effect of ∇B
 - With fixed plasma length and $B_{\min}$
5. Summary




SuSI - SUPERCONDUCTING SOURCE FOR IONS

- Fully superconducting with 6 solenoids
- Two 18 GHz klystrons, up to 3 kW
- Magnetic field up to 24 GHz
- Movable injection baffle, bias disk and puller electrode
- Extraction HV up to 27 kV
- Collimation channel, on-axis axial sputtering, on/off-axis axial oven




- More information: MOCOAK02 (L. T. Sun), MOCOBK01 (G. Machicoane), MOPOT16 (D. Cole)



X-RAY MEASUREMENT SETUP



X-RAY MEASUREMENT SETUP





X-RAY MEASUREMENT SETUP



- 2 collimators: Pb (83x83x83 mm) with W inserts (50 mm, OD 12 mm, hole $\varnothing$ 0.8 mm)
- Additional shielding:
 - Pb blocks inside the dipole
 - Stand-alone Pb shield between the collimators, hole $\varnothing$ 5 mm
 - ECR room Pb wall, hole $\varnothing$ 5 mm



X-RAY MEASUREMENT SETUP




- Ge detector located behind the wall
- Pb shielding around the detector to decrease background
- Collimators aligned through injection with laser
- 30 min measurement time
- Distance between extraction and detector: ≈ 3.2 m

DIFFERENT COLLIMATION GEOMETRIES

- Area seen by the detector at the extraction aperture ($\varnothing$ 10 mm) of SuSI
 - #1 $\varnothing$ 7 mm: All collimators
 - #2 $\varnothing$ 62 mm: 1st collimator removed
 - #3 $\varnothing$ 83 mm: #2 + Pb blocks inside the dipole & stand-alone Pb shield removed
 - #4 $\approx$ The whole ECRIS: #3 + 2nd collimator removed




DIFFERENT COLLIMATION GEOMETRIES



- Xenon/oxygen plasma, 1200 W
- Endpoint energies:
#1: ≈ 400 keV, #2: ≈ 650 keV, #3: ≈ 850 keV, #4: ≈ 1200 keV


DIFFERENT COLLIMATION GEOMETRIES



Difference in the shape:


- Plasma (cyclotron/synchrotron radiation) & collisions
- High energy electron losses
- Scattered radiation

EFFICIENCY CORRECTION AND PILE-UP REJECTION



- Energy spectra with and without efficiency correction and pile-up rejection for #1, #2, #3 and #4


EFFICIENCY CORRECTION AND PILE-UP REJECTION



1st figure: spectral temperature (fit range 130–175 keV)

no eff. corr. no pile-up reject.	eff. corr. no pile-up reject.	no eff. corr. pile-up reject.	eff. corr. pile-up reject.
56.8 ± 2.8 keV	49.4 ± 2.2 keV	56.6 ± 5.4 keV	49.2 ± 4.1 keV


EFFICIENCY CORRECTION AND PILE-UP REJECTION



1st figure: endpoint energies

no eff. corr. no pile-up reject.	eff. corr. no pile-up reject.	no eff. corr. pile-up reject.	eff. corr. pile-up reject.
425 ± 20 keV	440 ± 20 keV	385 ± 20 keV	375 ± 20 keV


EFFICIENCY CORRECTION AND PILE-UP REJECTION



4th figure: spectral temperature (fit range 600–800 keV)

no eff. corr. no pile-up reject.	eff. corr. no pile-up reject.	no eff. corr. pile-up reject.	eff. corr. pile-up reject.
118.9±0.6 keV	129.9±0.7 keV	117.7±1.3 keV	128.5±1.5 keV

EFFICIENCY CORRECTION AND PILE-UP REJECTION




4th figure: endpoint energies

no eff. corr. no pile-up reject.	eff. corr. no pile-up reject.	no eff. corr. pile-up reject.	eff. corr. pile-up reject.
1340±30 keV	1430±30 keV	1140±30 keV	1230±30 keV


FLATB VS. BMIN OPERATION

- FlatB operation
 - Unstable plasma
 - Very sensitive to tuning
 - Very high radiation levels at the extraction side (outside Pb shielding): 50 mR/h (FlatB) vs. 6 mR/h (Bmin)
- 0.1 T difference in Binj
- Tuned to $^{129}\text{Xe}^{27+}$, oxygen as a support gas
- 1500 W (18.0 GHz) + 350 W (17.8 GHz)





#profile	gradB (inj, ext)	plasma length	res. location (inj, ext)	$B_{\min}$
Bmin	-5.73 T/m, 5.01 T/m	124.5 mm	237.3 mm, 361.8 mm	0.47 T
FlatB	-2.19 T/m, 2.25 T/m	127.1 mm	229.1 mm, 356.1 mm	0.59 T

FLATB VS. BMIN OPERATION



- SuSI tuned in both cases keeping RF power constant
- $170 \mu\text{A}$ of Xe^{27+} with FlatB and Bmin
- CSD shifts towards lower charge states with FlatB
 - Reasons: injection field(?), $B_{\text{min}}(?)$, $\nabla B/\text{EEDF}(?)$, plasma-wave coupling(?)
 - FlatB: 4.4 mA, ext. pressure $8.2\text{e-}8$ torr
 - Bmin: 4.5 mA, ext. pressure $9.5\text{e-}8$ torr

FLATB VS. BMIN OPERATION




- Fit range for spectral temperature values: 130–200 keV


Spectral temperature		Endpoint energy	
FlatB	Bmin	FlatB	Bmin
68.1 ± 1.1 keV	58.0 ± 1.7 keV	850 ± 30 keV	450 ± 30 keV

∇B WITH FIXED PLASMA LENGTH AND $B_{\min}$

- 6 coils of SuSI allow tuning the ∇B while keeping profile almost unchanged
- Differences:
 - ∇B : 23 % (inj), 13 % (ext)
 - Plasma length: 1.1 mm
 - $B_{\min}$: 6.5 %
- Source tuned with profile 'med gradB' and then only the field profile is changed





#profile	gradB (inj, ext)	plasma length	res. location (inj, ext)	$B_{\min}$
med gradB	-5.73 T/m, 5.01 T/m	124.5 mm	237.3 mm, 361.8 mm	0.47 T
low gradB	-5.19 T/m, 4.75 T/m	125.4 mm	237.8 mm, 363.2 mm	0.49 T
high gradB	-6.38 T/m, 5.39 T/m	124.3 mm	236.1 mm, 360.4 mm	0.46 T



- SuSI tuned to Xe^{20+} with constant RF power (1500 W)
- CSD shifts to higher charge states while lowering ∇B
 - Reasons: $\nabla B/\text{EEDF}(\?)$, plasma-wave coupling(?)
 - low ∇B : 5.0 mA, ext. pressure $1.4\text{e-}7$ torr
 - med ∇B : 5.0 mA, ext. pressure $1.6\text{e-}7$ torr

∇B WITH FIXED PLASMA LENGTH AND $B_{\min}$



- Fit range for spectral temperature values: 130–200 keV

Spectral temperature

low ∇B	med ∇B	high ∇B
60.1 ± 0.9 keV	57.8 ± 0.8 keV	56.4 ± 0.9 keV


Endpoint energy

low ∇B	med ∇B	high ∇B
470 ± 20 keV	435 ± 20 keV	405 ± 20 keV

SUMMARY

When accurate spectroscopy information is needed:

- Efficiency correction...
- Pile-up rejection...
- Unfolding of the spectra (LBNL group)
 - Has impact on spectral temperature but not endpoint energy
- Energy resolution also affects the determination of endpoint energy (e.g. Ge vs. NaI)



SUMMARY

When accurate spectroscopy information is needed:

- Efficiency correction...
- Pile-up rejection...
- Unfolding of the spectra (LBNL group)
 - Has impact on spectral temperature but not endpoint energy
- Energy resolution also affects the determination of endpoint energy (e.g. Ge vs. NaI)


Different collimation geometries:

- By using different solid angles the shape and the endpoint energies of the X-ray spectra change significantly
- One can choose to study only the plasma or also higher energy electrons
- Shielding from the scattered radiation is very important

SUMMARY

FlatB operation:

- Same amount of medium charge states of xenon than with Bmin operation, but
- Very high radiation levels
- Endpoint energies almost doubled
 - possible issues with LHe consumption and radiation damage to insulator materials
- *However:* FlatB needs to be studied further...



SUMMARY

FlatB operation:

- Same amount of medium charge states of xenon than with Bmin operation, but
- Very high radiation levels
- Endpoint energies almost doubled
→ possible issues with LHe consumption and radiation damage to insulator materials
- *However:* FlatB needs to be studied further...

Changing ∇B with fixed field profile

- A clear impact on the heavy ion beam intensities
 - Smaller difference with oxygen
- A well-known trend in the X-ray spectra
- The effect of just the gradB shown now for the first time experimentally

ACKNOWLEDGMENTS

G. Machicoane (ion sources)
L. T. Sun (ion sources)
D. Cole (ion sources)
L. Tobos (ion sources)
A. Stolz (nuclear spectroscopy)
R. Fontus (mechanical design)

